

PRODUCT INFORMATION

Model: RLVF08204 and RLVF08204/LK (Lockable)

Type: Labcold Sparkfree Laboratory Freezer (-20°C)

For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future.



Image for illustrative purposes only

Capacity: 245 Litres (gross)

Construction: White enamel exterior, plastic interior.

Internal Arrangement: 6 drawers – wire mesh with transparent, cold retaining fronts.

Exterior Dimensions (mm): 1425 high x 595wide x 620 deep

Features:

Sparkfree Interior

CFC and HCFC free

Temperature set-point range: -18°C to -26°C

Analogue Temperature Display

Visual alarm for high temperature.

Top mounted ice packs to improve hold-over performance.

Low energy operation

Operator controlled defrost

Order Reference: RLVF 08204

RLVF 08204/LK