

## **PRODUCT INFORMATION**

**Model: RLVL03202**

**Type: Labcold Economy Laboratory Freezer**

For laboratory and general use where reasonable performance and low cost are important.



**Capacity: 90 Litres (gross)**

**Construction: White enamel exterior, fully moulded plastic interior.**

**Internal arrangement: Four drawers separated by fixed refrigerated shelves.**

**Exterior Dimensions (mm): 855 high x 505 wide x 500 deep**

### **Features:**

**Removable freezer drawers**

**CFC and HCFC free**

**Temperature set-point range: -15 to -20°C**

**Operator controlled defrost**

**Order Reference: RLVL 03202**