

PRODUCT INFORMATION

Models: RLFM0040, RLFM0075, RLFM0125, RLFM0165

Type: Labcold Flake Ice Maker

For the production and storage of flake (crushed) ice from tap water for laboratory purposes.



Construction: Stainless steel exterior, fully moulded plastic collection bin.

Service requirements: Mains water supply min 1 bar pressure at all times.

(Please bear in mind other demands on same water supply circuit)

240V at 50Hz, 3amp mains electricity

Drain of max height 150mm above floor level (if legs fitted).

Features:	RLFM0040	RLFM0075	RLFM0125	RLFM0165
Ice Production kg/24hrs*	40	75	125	165
Bin Storage Capacity kg	10	25	42	60
Height (with legs), mm	635	845**	1030**	915**
Width, mm	345	485	680	955
Depth, mm	520	570	588	590

Order Reference: RLFM 0040 RLFM 0075 RLFM 0125 RLFM 0165

*with inlet water and surrounding air temperatures both at 10°C

**can be reduced by 130 mm if supplied legs not fitted, but drain must then be at ground level.