

PRODUCT INFORMATION

Model: RLDF13043

Type: Labcold Pharmacy/Vaccine Refrigerator

For laboratory, surgery and general use for the storage of drugs and vaccines. Conforms in full with New UK Department of Health "Green Book" requirements.



Capacity: 360 litres (gross)

Construction: White enamel exterior, fully moulded plastic interior

Shelves Supplied: 6, Plastic coated wire mesh

Exterior Dimensions (mm): 1825 high x 600 wide x 615 deep

Features:

Digital display of temperature with max/min record

Audible alarm with mute facility

Fan assisted circulation

Door Lock

Easy-grip door handle

CFC, HCFC and Ammonia free

Accurate temperature control: +2°C to +8°C

Automatic Defrost

Order Reference: RLDF 13043