

PRODUCT INFORMATION

Model: LPTU0004

Type: Labcold Cytotherm dry-system rapid plasma thawer

For the rapid thawing of up to four packs (dependent upon size/pack type) of Fresh Frozen Plasma without direct water contact or the possibility of overheating.



Capacity: Up to 4 packs

Speed: Thaws from -35°C in approx 16 minutes

Construction: Moulded plastic interior and exterior

Exterior dimensions (mm): 330H x 332W x 520D

Features:

Indirect, circulating water system.

Digital display of time and temperature

Over-temperature alarm

Packs massaged by pulsing action

Replaceable bladders

Plasma packs remain separated from circulating water and dry at all times

Easy to clean and sterilise

Max temperature of water: 38°C

Order reference: LPTU 0004