

PRODUCT INFORMATION

Model: RLCG0600

Type: Labcold Advanced Refrigerated (cooled) Incubator (Glass door)

For maintaining samples at temperatures accurately controlled within the range +10 to 40°C. Temperatures and/or internal lighting can be cycled where appropriate accessories are fitted.



Capacity: 556 litres (gross)

Construction: White enamel exterior, stainless steel interior. Mounted on castors. Glass door

Shelves supplied: 4; Slatted stainless steel

Exterior dimensions (mm): 1995 high x 730 wide x 800 deep

Features:

Microprocessor temperature control with digital display

Temperature set-point range: +10 to 40°C

Temperature precision: better than +/- 0.5°C at 20°C

Fan circulated interior

Door lock

CFC and HCFC free

Manual defrost with condensate drain

Order reference: RLCG 0600