

## PRODUCT INFORMATION

**Model: RLCH0600**

**Type: Labcold Advanced Refrigerated (cooled) Incubator**

For maintaining samples at temperatures accurately controlled within the range 0 to 40°C\*. Temperatures and/or internal lighting can be cycled where appropriate accessories are fitted.



**Capacity: 556 litres (gross)**

**Construction: White enamel exterior, stainless steel interior. Cabinet mounted on castors.**

**Shelves supplied: 4; Slatted stainless steel**

**Exterior dimensions (mm): 1995 high x 730 wide x 800 deep**

### Features:

**Microprocessor temperature control with digital display**

**Temperature set-point range: 0 to 40°C\***

**Temperature precision: better than +/- 0.5°C at 5°C**

**Fan circulated interior**

**Door lock**

**CFC and HCFC free**

**Manual defrost with condensate drain**

**Maximum ambient temperature: 32°C**

**\*Alternative temperature ranges available to special order**

**Order reference: RLCH 0600**