

SPARKFREE CHEST FREEZERS



RLCF0720



RLCF2120

Chest models are a viable alternative, where space permits, as the temperature rise when the door is opened is less than upright counterparts. Locating samples can be made simple and straightforward by use of a Labcold sample racking system (see page 21).

Part Number	RLCF0720	RLCF1520	RLCF2120
Configuration	Chest	Chest	Chest
Capacity (litres)	215	447	607
Exterior size (mm)			
High	860	860	860
Wide	720	1300	1700
Deep	730	730	730
Sparkfree interior	✓	✓	✓
Internal fan	✗	✗	✗
Door lock	✓	✓	✓
Temp display	✗	✗	✗
Alarm	✓	✓	✓
Auto defrost	✗	✗	✗
Low energy	✓	✓	✓
No. shelves	N/A	N/A	N/A
Operating temp (nominal)	-20°C	-20°C	-20°C

✓ – standard ○ – available as an optional accessory ✗ – not available on this model

**ILSHIN
ULTRA-LOW
TEMPERATURE
FREEZERS**

Page 18-19

