

PRODUCT INFORMATION

Model: RLHD13043 and RLHD13043/LK (lockable)

Type: Labcold Sparkfree Laboratory Refrigerator

For laboratory and general use where there is a concern that volatile, flammable materials or components might be stored now, or at some point in the future.



Capacity: 360 litres (gross)

Construction: White coated steel, fully moulded plastic interior

Shelves Supplied: 4; Plastic coated wire mesh

Exterior Dimensions (mm): 1755 high x 585 wide x 615 deep

Features:

Sparkfree Interior

Digital temperature display

Easy-grip door handle

CFC and HFC free

Temperature set-point range: +2 to +8°C

Automatic Defrost

Please note: In order to provide a spark free interior no internal air-circulating fan is fitted to this product. Some gradient in temperature between top and bottom shelves is therefore to be expected.

Order Reference: RLHD 13043

RLHD 13043/LK