

# SPARKFREE LABORATORY FREEZERS



RLVF04202



RLVF1220



RLVF08204

- Sparkfree safety interiors
- Lockable models
- CFC, HCFC and Ammonia free
- Operating temperature  $-20^{\circ}\text{C}$
- Low energy models
- Operator controlled defrost
- Maximum ambient temperatures up to  $32^{\circ}\text{C}$

Laboratory safety freezers, with sparkfree interiors, to match the Labcold sparkfree or Labcold fan circulated refrigerator ranges. A choice of front opening and chest configurations is offered with all models featuring operator controlled defrost for improved temperature stability. Locks, alarms and other accessories can be fitted at time of manufacture to most models where not featured as standard.

Options available (dependant on model): models with door lock; temperature alarm/display; wall mounting kit for 45 litre models; sample storage racking.

For alarm/display details see page 27.

Part Number	RLVF 02202	RLVF 02202/LK	RLVF 04202	RLVF 04202/LK	RLVF 08204	RLVF 08204/LK	RLVF 10202	RLVF 10202/LK	RLVF 1220	RLVF 1220/LK	RLVF 14201
Configuration	Counter top	Counter top lockable	Under-bench	Under-bench, lockable	Free-standing	Free-standing lockable	Free-standing	Free-standing lockable	Free-standing	Free-standing lockable	Free-standing
Capacity (litres)	45	45	120	120	245	245	296	296	344	344	375
Exterior size (mm)											
High	525	535	850	850	1435	1435	1750	1750	1850	1850	1850
Wide	500	500	550	550	600	600	604	604	595	595	595
Deep	500	500	620	620	600	600	670	670	595	595	640
Sparkfree interior	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Internal fan	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗
Door lock	✗	✓	✗	✓	✗	✓	✗	✓	✗	✓	✓
Temp display	○	○	○	○	○	○	✓	✓	✓	✓	✓
Alarm	○	○	○	○	○	○	○	○	○	○	○
Auto defrost	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗	✗
Low energy	✗	✗	✓	✓	✗	✗	✓	✓	✗	✗	✗
No. shelves/drawers	1	1	4	4	6	6	9	9	5/3	5/3	6
Operating temp (nominal)	$-20^{\circ}\text{C}$	$-20^{\circ}\text{C}$	$-20^{\circ}\text{C}$	$-20^{\circ}\text{C}$	$-20^{\circ}\text{C}$	$-20^{\circ}\text{C}$	$-20^{\circ}\text{C}$	$-20^{\circ}\text{C}$	$-20^{\circ}\text{C}$	$-20^{\circ}\text{C}$	$-20^{\circ}\text{C}$

✓ - standard    ○ - available as an optional accessory    ✗ - not available on this model