

SUPERFREEZERS



RLBF1740



RLHE1145



RLHD14301

- **-30°C to -45°C units**
- **Microprocessor Control**
- **Audible and visual temperature alarm**
- **Sparkfree and fan assisted models**
- **Upright or chest models**

Superfreezers is the term Labcold uses to describe freezers designed to operate at temperatures below the more common -20°C but are not as cold as Ultra Low Temperature -80°C freezers.

Front opening models are offered in bench top -40°C, and free-standing -30°C configurations, a further three sizes of top opening (chest) -45°C freezer are also available.

Options available: remote alarm contacts; telephone dial-out unit; sample storage racking.

Part Number	RLBF1740	RLBF1740/LK	RLHD14301	RLHE0845	RLHE1145	RLHE1345
Configuration	Benchtop	Benchtop, lockable	Free-standing	Chest	Chest	Chest
Capacity (litres)	50	50	400	237	314	375
Exterior size (mm)						
High	580	580	1995	860	860	860
Wide	550	575	600	1050	1305	1500
Deep	620	620	700	735	735	735
Sparkfree interior	✓	✓	✗	✓	✓	✓
Internal fan	✗	✗	✓	✗	✗	✗
Door lock	✗	✓	✓	✓	✓	✓
Temp display	✓	✓	✓	✓	✓	✓
Alarm	✓	✓	✓	✓	✓	✓
Auto defrost	✗	✗	✗	✗	✗	✗
No. shelves/baskets	1	1	4	1	2	2
Operating temp (nominal)	-40°C	-40°C	-30°C	-45°C	-45°C	-45°C

✓ - standard ✗ - not available on this model