

ORDER BASIC MODEL L46 TOGETHER WITH ADDITIONAL ACCESSORIES BELOW



A

 B
b1

 C

 D

 E

A COMPLETE MIXER-LINE FROM ONE BASIC UNIT

A	MICRO-TITRE-PLATE -HOLDER	easy adaptable holder for one 96-well micro-titre-plate shaking with Model L 46
	Ordering number	46.100
B	MINI-REACTION-TUBE-HOLDER	easy adaptable holder for 20 mini-reaction-tubes; including polythene-foam holder
	Ordering number	46.200
	b1 Spare polythene-foam-holder for 20 mini reaction-tubes	
	Ordering number	46.250
C	One-position TEST-TUBE-HOLDER	easy adaptable holder for a single test-tube (max. 20x200 mm) Including aluminium rod(M10)
	Ordering number	46.300
D	Multi-position TEST-TUBE-HOLDER	easy adaptable holder for up to 7 test-tubes (18x115 mm) (max) Changing while mixing is possible
	Ordering number	46.400
E	Universal CONTACT-PLATE	easy adaptable universal contact-plate for touching and mixing up to 6 tubes at once.
	Ordering number	46.500