

MODEL LD-79

TEST-TUBE-ROTATOR



50/60 HZ CAT.NO. 79000 (230V) / 79010 (115V)

ELECTRONICALLY CONTROLLED STEPPER-MOTOR PROVIDES:

- Digital read-out for - speed 5-50 rpm.
- time 0-99 hrs/1-59 min.
- Electronic timer for "auto-stop".
- "Memory-function" for repeated operation.
- "Re-start-function" for continued operation after line-voltage failure including auto-stop.
- "Settings saving" available.



50/60 HZ CAT.NO. 28000 (230V) / 28010 (115V)

MODEL L-28 TEST-TUBE ROTATOR

- Provided with strong 10 rpm. synchronic motor
- Accepts all disks mentioned in the listing besides
- Dimensions (LxWxH) 210x200x220 mm.
- Weight: 3,5 Kg.

- Provides mixing for diffusion, dialysing, extraction etc.
- Tube diameter range 4-40 mm.
- Accepts all disks mentioned in the listing below.
- Mixing action can be set by:
 - speed control 5-50 rpm.
 - angle of rotation 90° to 27°.
- Quick/easy exchange of disks by unscrewing central nut.
- Dual or triple disk support may be carried out.
- Tubes may be exchanged while the rotator is running.
- Dimensions: 200 (L) x 200 (W) x 250 (H) mm.
- Weight: 3,5 Kg.



ROTATING DISKS				FOR TEST TUBES		
FOR MINI-REACTION VIALS/TUBES				CATALOG NUMBER	NUMBER CLAWS	TUBE Ø RANGE
CATALOG NUMBER	NUMBER CLAWS	TUBE Ø RANGE	VOLUME ML			
28006	40	4-6	0.4	28020	20	12-16
28009	30	6-9	0.5	28030	20	16-19
28012	20	9-12	1.5	28032	10	28-32
	28090	adapter for dual disk support		28040	10	32-40
				28050	15	22-25
				28070	17	19-22
				28080	12	25-28



MODEL LD-76

ROTARY MIXER

- Provides serial rotation of tubes and mini-tubes for extraction, diffusion, etc.
- Mixing action may be set by speed control 5-50 rpm.
- Tube-diameter-range 4-32 mm. Chosen from listing below.
- Quick/easy coupling/exchange of rotation drums/bars.
- Tubes may be exchanged while the rotator is running.
- Dimensions: LxWxH 450x170x190 mm.
- Weight: 3,3 Kg.

50/60 HZ CAT.NO. 76000 (230V) / 76010 (115V)

ELECTRONICALLY CONTROLLED STEPPER-MOTOR PROVIDES:

- Digital read-out for : - speed 5-50 rpm.
- time 0-99 hrs/1-59 min.
- Electronic timer for "auto-stop".
- "Memory-function" for repeated operation.
- "Re-start-function" for continued operation after line-voltage failure including auto-stop.
- "Settings saving" available.


 ON/OFF MAINS

 START/STOP MOTOR


 "MODE" SWITCH

 "UP/DOWN"





 "TIMER/CYCLING"



LABINCO



ACCESSOIRES FOR ROTARY-MIXERS

Bar/Drum outfit	tube \varnothing range mm	cat.no. for L26/LD-76	max.tube lenght mm	max no. tubes
	16-20	26020	190	20
	19-22	26022	185	18
	22-27	26027	180	18
	4-6	26036	85	104
	6-9	26039	85	104
	9-12	26040	85	104
	12-16	26480	120	48
	if 2 rows used		140	24
	9-12	26560	120	56
	if 2 rows used		140	28