

211DS Shaking Incubator

Designed for a wide variety of applications



- Exclusive SmartChek temperature control system
- Stackable design saves space
- Full range of accessories available
- Internal electrical outlet

Product Information

A small footprint, compact, stackable design and expanded temperature range make the 211DS incubators ideal for molecular biology and general use. Its integral orbital shaker makes it highly versatile.

Labnet's exclusive SmartChek temperature control system guarantees precise temperature control. Mechanical convection maintains a stable temperature environment and provides quick recovery after opening the door. A safety thermostat is located on the back of the unit. The insulated door has a large, double glass observation panel and opens nearly 180° for easy access.

Two full and one half shelf are supplied with each incubator. A flat platform and predrilled flask platform (clamps sold separately) are available for the 211DS

Internal electrical outlet allows using small equipment such as GyroMini or Mini LabRoller Rotator inside the chamber of the 211DS.



Bottles



Mini LabRoller

Product Specifications

Chamber volume	1.7 cu. ft.
Temperature range/display	30° to 95°C ¹ in 0.1° increments/digital
Temp. uniformity/accuracy	±0.25°C/±0.1°C
Shaker speed/orbit (DS,)	Variable, 20 to 400rpm/19 mm
Flask capacity (DS,)	4 x 1 L, 5 x 500 ml, 9 x 250 ml, 16 x 125 ml
Construction	Interior - mirrored stainless, exterior - cold rolled steel
Chamber (W x D x H)	13.5 x 14.75 x 15 inches/34.3 x 37.5 x 38.1 cm
Exterior (W x D x H)	16.75 x 21.75 x 23 inches/42.5 x 55 x 58 cm
Weight	D -73 lb/32.2 kg, DH- 83 lb/37.7 kg, DHS - 85 lb/38.2 kg
Electrical*	120V~, 60 Hz

¹Maximum 80°C in 211DS *Other voltages available upon request

Ordering Information

Catalog Number	Item Description
I-5211-DS	Digital Incubator/Shaker, Model 211DS
I-5211-SA	Stacking adapter, required when stacking incubators
I-5230	Flask clamp platform, predrilled, fits 211DS, 211DHS
I-5231	Flat platform with nonslip rubber mat, fits 211DS, 211DHS
I-5231-D	Double flat platform with nonslip rubber mat, fits 211DS, 211DHS
S2040-85	Flask clamp, 125 ml (max. 15)
S2040-99	Flask clamp, 250 ml (max. 9)
S2040-09	Flask clamp, 500 ml (max. 5)
S2040-01	Flask clamp, 1000 ml (max. 4)