

Wolf's Drying Cabinet Range



Features

- Maximum temperature 50°C above ambient
- Energy regulator control
- Easy clean stove powder paint finish
- Toughened glass doors
- Chrome plated steel wire shelves
- 113 & 226 litre fitted with sliding doors
- 540 & 1000 litre fitted with hinged doors
- 540 & 1000 litre available with fan circulation and or variable rate extraction

Optional Extras

- Maximum temperature of 65°C above ambient is also available on the 540 & 1000 litre Drying Cabinets
- Stainless steel body available on 113 & 226 litre units

Design

A range of cabinets for drying drained glassware and instruments, or for applications requiring controlled warming. The construction of the cabinets are of single skin mild steel and finished in easy to clean stove powder paint for durability. The chromium plated grid wire shelves are adjustable on a ladder racking system. The cabinets are base heated by incoloy sheathed elements and are powered by 220/240v single-phase supply.

DC series	113	226	540	540/F	540/EX
Ext. Dims (mm) H x W x D	559 x 730 x 360	666 x 928 x 449	1664 x 583 x 610	1814 x 583 x 610	1924 x 583 x 685
Int. Dims (mm) H x W x D	460 x 728 x 341	567 x 926 x 430	1550 x 581 x 600	1550 x 581 x 600	1550 x 581 x 600
Volume Approx. (litres)	113	226	540	540	540
Number of Shelves	3	3	3	3	3
Shelf Size (mm) W x D	810 x 330	910 x 415	565 x 585	565 x 585	565 x 585
Power Rating (Watts)	750	1000	1750	1750	1750
Net Weight (Kg)	30	40	80	83	84

DC series	540/F/EX	1000	1000/F	1000/EX	1000/F/EX
Ext. Dims (mm) H x W x D	1924 x 583 x 685	1664 x 1110 x 585	1814 x 1110 x 585	1924 x 1110 x 660	1924 x 1110 x 660
Int. Dims (mm) H x W x D	1550 x 581 x 600	1550 x 1108 x 575	1550 x 1108 x 575	1550 x 1108 x 575	1550 x 1108 x 575
Volume Approx. (litres)	540	1000	1000	1000	1000
Number of Shelves	3	4	4	4	4
Shelf Size (mm) W x D	565 x 585	1095 x 665	1095 x 665	1095 x 665	1095 x 665
Power Rating (Watts)	1750	2750	2750	2750	2750
Net Weight (Kg)	87	110	114	114	118