

## Undercounter Laboratory Fridge & Freezer

Spark-free and lockable, specifically designed for under counter laboratory use.



**830mm  
high**

**Lec+**

## Undercounter Laboratory Fridge & Freezer

Designed specifically for use in laboratories, this range is just 830mm high so will fit neatly underneath the work bench. Being spark-free and lockable, the units offer safe and reliable storage of volatile liquids, reagents, chemicals and samples.

### Control and Efficiency

- Easy to use integral control panel makes checking the internal temperature straightforward.
- Internal air temperature is kept stable as there is no need to open the door to check the temperature.
- Automatic defrost on the LR307 ensures more efficient cooling.

### Security

- Spark-free interior reduces the risk of an internal explosion.
- Lockable door to ensure only authorised individuals are able to access the contents.
- An audio-visual temperature alarm sounds if the internal temperature goes outside the set alarm limits, helping reduce the risk of ruining the contents of the fridge/freezer.

### Flexibility

- Reversible doors give greater flexibility when positioning the product.
- Removable shelves on the LR307 give maximum storage capacity.

### Environmentally Friendly

- Low energy consumption.
- CFC and HCFC free.

Model	LR307
SKU number	444449326
Product dimensions (in mm)	495w x 580d x 830h
Packed dimensions (in mm)	560w x 655d x 850h
Product weight (in Kg)	31.5
Packed weight (in Kg)	35
Net capacity (in litres)	103
Gross capacity (in litres)	107
Temperature range	+2°C to +10°C
Colour	White
Door type	Hinged, reversible
Door style	Solid
Lock	Yes
Adjustable feet	Front Yes (1) Back Fixed
Temperature display in the panel	Yes
Temperature alarm	Yes
Shelves	Yes (2)
Wire drawer	Yes (1)
Interior light	No
Fridge defrost	Automatic
Fan	No
Control	Electronic
Climate class	N
Coolant	R134a
AUX Output	No
UK 3-pin plug	Yes
Voltage supply	230V / 50Hz

Model	ISU37
SKU number	444449327
Product dimensions (in mm)	495w x 580d x 830h
Packed dimensions (in mm)	560w x 655d x 850h
Product weight (in Kg)	33
Packed weight (in Kg)	36.5
Net capacity (in litres)	87
Gross capacity (in litres)	102
Temperature range	-18°C to -20°C
Colour	White
Door type	Hinged, reversible
Door style	Solid
Lock	Yes
Adjustable feet	Front Yes (1) Back Fixed
Temperature display in the panel	Yes
Temperature alarm	Yes
Shelves	Yes (2)
Wire drawer	Yes (1)
Interior light	No
Freezer defrost	Manual
Fan	No
Control	Electronic
Climate class	N
Coolant	R134a
AUX Output	No
UK 3-pin plug	Yes
Voltage supply	230V / 50Hz

All dimensions are approximate.  
We reserve the right to change the specifications of any product without prior notice.

**Lec**+

Working Together