

AXL 230-240V/50-60 Hz



Leica Autostainer XL (ST5010)-- Automated Slide Stainer

045635136

Leica Autostainer XL (ST5010) - Automated Slide Stainer
230-240 V / 50-60 Hz

The instrument features:

Fully self-contained bench top unit.

18 reagent vessels, 1 integrated fan-forced oven, 5 wash vessels with flow control valves, 10 slide racks each holding 30 standard slides. Wide range of rack adapters.

Continuous load via load and exit drawers. Menu-driven, user-friendly integrated control panel.

15 programs 25 steps each, ability to run multiple programs together, station sequence, incubation time and agitation flexible.

Air filter vent with charcoal filter to retain solvent fumes.

Communication and remote alarm port.

Designed and manufactured in compliance with UL/c-UL and CSA regulations.

Standard delivery includes:

1 Basic instrument (230-240 V/50-60 Hz)

1 Accessory-Kit (0456 35660) consisting off:

-22 Reagent vessels with lids (0475 33659)

-5 Wash vessels (0456 35268)

-10 Slide racks 30, plastic (0475 33750)

-2 Slotted lids for reagent vessels (0475 34486)

-1 Jumper cable, mains (0411 34604)

-1 1 Tubing band clamp (0422 31972)

-1 Elbow hose (0475 33669)

-1 V-filter 3/4" for inlet hose (0456 36101)

-1 1 Jack connector (remote alarm) (6844 01005)

1 Activated carbon filter (0474 32273)

1 Fume cover, front (0456 35240)

1 Wax collection tray (heater)

(included in the instrument) (0456 35216)

1 Program pad (attached to the instrument) (0456 35459)

1 Inlet hose with washer (0474 32325)

1 Drain hose (0475 35748)

1 Set of power cords:

-1 Power cord "UK" ST-BU F5-A (0411 27822)

-1 Power cord "D" (0411 13558)

1 Instruction manual Leica ST5010 - G/E/F/S including 1 reference manual

