

Leica EG1150 C Cold Plate



230 V/50 Hz Part No 038838037

Cold plate for approx. 70 molds with easy-to-clean surface.
Illuminated ON / OFF switch on the front side of the system.
Controlled, constant temperature of -5°C.

Technical Data:

General instrument specifications:

- Admissions: CE
- Mains voltage 230 V/50 Hz
- Maximum power draw: 350 VA
- Maximum start-up current (5s): 25 A
- Mains fuse 230 V:- overload cutout T 4 A
- Operating temperature range:+18°C to +40°C
- Relative air humidity: max. 60%, non-condensing

Cooling system:

- Cooling capacity: 110 W
- Refrigerant: R134a

Dimensions and weight:

- Capacity: approx. 70 molds
- Size of Cold Plate (D x W): 370 x 330 mm
- Size of Cold Plate (D x W): 14.57 x 12.99 inches
- Size of Instrument:
- Width: 360 mm / 14.17 inches
- Depth: 650 mm / 25.59 inches
- Height: 360 mm / 14.17 inches
- Weight: 23 kg / 50.66 lbs

Standard delivery includes:

- 1 Basic instrument (0388 38037)
- 1 Set of Power cords:
- Power cord "D" (0411 13558)
- Power cord "UK"ST/BU F-5A (0411 27822)
- 1 Fuse 6.3 x 32mm T1A (6943 01001)
- 1 Instruction Manual Leica EG1150 C

--	--