

Low Capacity Chromatography Cabinets

These compact chromatography cabinets are manufactured by Liebherr.

One access ports are supplied – state left or right at point of order – to allow cables and tubes to be introduced into the cabinet. Temperature is maintained by means of a recirculating fan and toughened insulated glass door.

- Recirculated air cooling
- Automatic defrost
- Lockable glass door
- Left- or right-hand opening
- Digital thermometer
- Interior light
- Six shelves



Specifications		
Model	UKS3602M	UKS5002M
Gross capacity, L	352	491
Temperature range, °C	+2 to +15	+2 to +15
Internal dimensions (h x w x d), mm	1418 x 513 x 433	1338 x 646 x 516
External dimensions (h x w x d), mm	1589 x 602 x 600	1515 x 752 x 720
Weight, kg	77.5	93.0