

Liebherr GG5210-740 Sparkfree Laboratory Freezer



Technical Data

Gross capacity, total [l]	513
Dimensions [cm]	H/W/D: 172.5/75/76
Net depth [cm]	56
Net width of interior compartment [cm]	60.7
Net height [cm]	151
Energy consumption in 24h [kWh]	1,2
Cooling system, freezer compartment	static
Defrosting, freezing compartment	manual
Temperature range freezer compartment [°C]	-14°C to -28°C
Housing material	steel
Colour	white
Door/lid material	steel
Material of interior containers	Plastic white
Type of control	Electronic control system
Temperature display, freezer compartment	external digital

Malfunction: Warning signal	optical and acoustical
potential-free contact	-
Storage shelves, freezer compartment	7
Storage shelf material, freezer compartment	-
Load capacity storage shelves, freezer compartment [kg]	60
Handle	Slimline door handle with integrated opening mechanism
Pedal opener	-
Height of adjustable feet [cm]	-
lock	-
Self-closing door	-
Door hinges	right interchangeable
Gross weight [kg]	111
Net weight [kg]	105
Climate class	SN-ST
Refrigerant	R 600a
connection rating [W]	1,0 A / 150
Voltage [V]	220-240