

Refrigerators with spark-free interiors

FKEX 1800, 2600, 3600 and 5000



These 4 models have been independently tested and certified in accordance with EU guideline 94/9/EG - ATEX 3500 X - effective from 1. July 2003



FKEX 1800

FKEX 2600

FKEX 3600

FKEX 5000

LIEBHERR

Die professionelle Kältetechnik.

FKEX models with spark-free interior.

Specifications:

- **Static cooling system**
- **Temperature range: +2°C to +10°C**
- **Adjustable glass shelves**
- **Plugged defrost water drain line**
- **Removable internal drip-tray**
- **External thermostat**
- **Warning notice in 2 languages (English/German) to door and inner liner indicating explosion-proof interior**
- **Operating instructions adapted to comply with the EU regulation**
- **One-piece moulded inner liner for easy hygiene**
- **Adjustable front legs**
- **Easy door re-hinge, lockable door**
- **Removable door gasket**
- **Suitable for ambient temperature from +10°C to +32°C**
- **Environmentally-friendly insulation and refrigeration system**



Technical data.

	Gross capacity (l)	Ext. dims. (H/W/D: mm)	Int. dims. (H/W/D: mm)	Voltage (V)	Energy consumption (kWh/24h)
FKEX 1800	180	885 / 600 / 600	734 / 513 / 441	230 / 50	0.90
FKEX 2600	260	1216 / 600 / 600	1047 / 513 / 441	230 / 50	1.00
FKEX 3600	360	1590 / 600 / 600	1418 / 513 / 441	230 / 50	1.10
FKEX 5000	500	1516 / 752 / 710	1338 / 652 / 537	230 / 50	1.10

LIEBHERR

Die professionelle Kältetechnik.