

FOCUS 1

Monocular Biological Microscope



Specifications	
Optical head	monocular, inclined at 45° and 360° rotatable
Eyepiece	10x/15mm wide-field with pointer, mounted student proof
Nosepiece	quadruple, revolving with click stop
Objectives	achromatic - 4x/0.10/160/- - 10x/0.25/160/0.17 - 40x/0.65/160/0.17 (spring loaded)
Precondenser with blue lens	
Stage	size 110x110mm with built-in condenser, iris diaphragm and filter holder. Two stage clamps. Stage predrilled to accept an attachable mechanical stage.
Focusing knobs	coaxial fine and coarse
Mains power	220-240V/50-60Hz (euro connector)
Converter	built-in low voltage (12V - electronic) with light intensity control
Bulb	halogen 12V/20W/G4
Supplied with	1 spare bulb, 2 spare fuses, dust cover and instruction manual
Packing	two styrofoam shelves in a cardboard box
Packing dimensions	28x19x37cm (WxDxH). Gross weight: 4 kg