

# COOL AND SAFE STORAGE


## Racks for chest freezers – stainless steel, sold by the each

Article no.	For boxes	HxDxW mm	No. of boxes
TE21450	50 mm	335x140x140	6
TE21452	50 mm	555x140x140	10
TE21454	50 mm	610x140x140	11
TE21456	50 mm	665x140x140	12
TE21478	50 mm	720x140x140	13
TE21458	75 mm	326x140x140	4
TE21460	75 mm	406x140x140	5
TE21462	75 mm	486x140x140	6
TE21464	75 mm	551x140x140	7
TE21466	75 mm	629x140x140	8
TE21465	85 mm	607x140x140	7
TE21468	100 mm	320x140x140	3
TE21470	100 mm	425x140x140	4
TE21472	100 mm	530x140x140	5
TE21474	100 mm	635x140x140	6



## Racks for upright freezers – stainless steel, sold by the each

Article no.	For boxes	HxDxW mm	Boxes (HxD)
TE21480	50 mm x 12	224x423x139	4x3
TE21485	50 mm x 16	224x562x139	4x4
TE21490	50 mm x 20	279x562x139	5x4
TE21575	50 mm x 24	331x562x139	6x4
TE21495	75 mm x 6	164x423x139	2x3
TE21500	75 mm x 8	164x562x139	2x4
TE21505	75 mm x 12	224x562x139	3x4
TE21510	100 mm x 6	214x423x139	2x3
TE21515	100 mm x 8	214x562x139	2x4



## Cryo-vials in polypropylene for temperatures down to -190°C

1.0 and 2.0 ml vials will fit into 50 mm boxes - 4.0 ml vials will fit into 75 mm boxes - 5.0 ml vials will fit into 100 mm boxes. All tubes are DNase and RNase free.

Order no.	Description	Bags of	Package
TE78005	1.0 ml, print, sterile	100	1000
TE78010	1.0 ml, with loose lids	100	1000
TE78020	2.0 ml, print, sterile	100	1000
TE78025	2.0 ml, with loose lids	100	1000
TE78035	4.0 ml, print, sterile	100	500
TE78040	4.0 ml, with loose lids	100	1000
TE78050	5.0 ml, print sterile	100	500
TE78055	5.0 ml, with loose lids	100	1000
TE78200	Rack 40 well for Cryo-vials		1



## Colour inserts for Cryo-vials

Order no.	Description	Package
TE78100	White	500
TE78101	Blue	500
TE78102	Yellow	500
TE78103	Pink	500
TE78104	Red	500
TE78105	Green	500
TE78106	Assorted	500

# COOL AND SAFE STORAGE

## Cardboard boxes WxD 133x133. Boxes and dividers are sold in packages of 36 each

Height	Article no. Box 50 mm	Article no. Box 75 mm	Article no. Box 100 mm
White	TE22005	TE22050	TE22100
Blue	TE22010	TE22055	TE22101
Red	TE22015	TE22060	TE22102
Green	TE22020	TE22065	TE22103
Yellow	TE22025	TE22070	TE22104

Grid dividers, height	Article no. 30 mm for 50 mm box	Article no. 40 mm for 75 mm box	Article no. 65 mm for 100 mm box
7x7 cells for vials up to 16 mm O.D.	TE22106	TE22130	TE22140
8x8 cells for vials up to 14 mm O.D.	TE22108	TE22132	TE22142
9x9 cells for vials up to 13 mm O.D.	TE22110	TE22134	TE22144
10x10 cells for vials up to 12 mm O.D.	TE22120	TE22136	TE22146



## Cardboard boxes WxD 136x136. Boxes and dividers are sold in packages of 36 each

Height	Article no. Box 50 mm	Article no. Box 75 mm	Article no. Box 100 mm
White	TE22500	TE22525	TE22540
Blue	TE22501	TE22526	TE22541
Red	TE22502	TE22527	TE22542
Green	TE22503	TE22528	TE22543
Yellow	TE22504	TE22529	TE22544

Grid dividers, height	Article no. 30 mm for 50 mm box	Article no. 40 mm for 75 mm box	Article no. 65 mm for 100 mm box
7x7 cells for vials up to 16 mm O.D.	TE22686	TE22706	TE22736
8x8 cells for vials up to 14 mm O.D.	TE22685	TE22705	TE22735
9x9 cells for vials up to 13 mm O.D.	TE22684	TE22704	TE22734
10x10 cells for vials up to 12 mm O.D.	TE22683	TE22703	TE22733




## Polypropylene boxes WxD 130x130 including grid divider 9x9. 50 and 80 mm boxes are sold in packages of 36 each. 90 mm boxes are sold in packages of 30 each

Height	Article no. Box 50 mm	Article no. Box 80 mm	Article no. Box 90 mm
Blue	TE22300	TE22325	TE22351
Green	TE22305	TE22330	TE22353
Clear	TE22310	TE22335	TE22350
Red	TE22315	TE22340	TE22352
Yellow	TE22320	TE22345	TE22354

## Polypropylene boxes WxDxH 130x130x130 mm – for 15 & 50 ml centrifuge tubes. Boxes are sold in packages of 12 each

Tubes (no. of tubes per box)	Article no. 15 ml (25)	Article no. 50 ml (10) + 15 ml (2)
Blue	TE22390	TE22360
Green	TE22391	TE22361
Clear	TE22392	TE22362
Red	TE22393	TE22363
Yellow	TE22394	TE22364



**TENAK**