

CM-Series (With Temp/ Controller)

Controlled Heating Mantles (analogue)

Specifications:

Case Construction: Steel, with chemical resistant coating
Thermal Insulation: Moulded silica powder
Heating Element: Nichrome wire
Mantle: Woven glass fibre mat
Element Temperature: Up to 400°C (Max)
Sensor Probe:: CA (K) Type, Ø 4.8mm * 200mm (1.5m lead)
Power Supply: 220V AC
Temp. Range (options): Low (up to 100°C) / High (up to 200°C)



MODEL NO.	CAPACITY	WATTS
CM-7005	50ml	60 W
CM-701	100ml	80 W
CM-702	250ml	150 W
CM-703	500ml	260 W
CM-704	1000ml	350 W
CM-705	2000ml	450 W
CM-706	3000ml	550 W
CM-707	5000ml	700 W
CM-708	10000ml	1200 W
CM-709	20000ml	1800 W