

C-Series (Without Temp/ Controller)

Un-Controlled Heating Mantles

Specifications:

Case Construction: Steel, with chemical resistant coating
 Thermal Insulation: Moulded silica powder
 Heating Element: Nichrome wire
 Mantle: Woven glass fibre mat
 Element Temperature: Up to 450°C (Max)
 Power Supply: 220V AC



MODEL NO.	CAPACITY	WATTS
C-4005	50ml	60 W
C-401	100ml	80 W
C-402	250ml	150 W
C-403	500ml	260 W
C-404	1000ml	350 W
C-405	2000ml	450 W
C-406	3000ml	550 W
C-407	5000ml	700 W
C-408	10000ml	1200 W
C-409	20000ml	1800 W