


DM-Series (With Temp/ Controller)

Controlled Heating Mantles (digital)

Specifications:

Case Construction: Steel, with chemical resistant coating
 Thermal Insulation: Moulded silica powder
 Heating Element: Nichrome wire
 Mantle: Woven glass fibre mat
 Element Temperature: Up to 400°C (Max)
 Sensor Probe:: CA (K) Type, Ø 4.8mm * 200mm (1.5m lead)
 Power Supply: 220V AC
 Temp. Range (options): Low (up to 100°C) / High (up to 200°C)



MODEL NO.	CAPACITY	WATTS
DM-6005	50ml	60 W
DM-601	100ml	80 W
DM-602	250ml	150 W
DM-603	500ml	260 W
DM-604	1000ml	350 W
DM-605	2000ml	450 W
DM-606	3000ml	550 W
DM-607	5000ml	700 W
DM-608	10000ml	1200 W
DM-609	20000ml	1800 W