

## EAM-Series (Extraction Mantle)

### Controlled Heating Mantles

#### Specifications:

Case Construction:	Body - Steel with chemical resistant coating Cover - High grade stainless steel
Thermal Insulation:	Moulded silica powder
Heating Element:	Nichrome wire
Mantle:	Woven glass fibre mat
Element Temperature:	Up to 450°C (Max)
Number Of Positions:	3 or 6 place models
Power Supply:	220V AC
Support Rod:	Ø 12.7mm



MODEL NO.	NO. POSITIONS	CAPACITY	WATTS
EAM-9201-03	3	100 ml	80 W (x3)
EAM-9202-03	3	250 ml	150 W (x3)
EAM-9203-03	3	500 ml	260 W (x3)
EAM-9204-03	3	1000 ml	350 W (x3)
EAM-9205-03	3	2000 ml	450 W (x3)
EAM-9201-06	6	100 ml	80 W (x6)
EAM-9202-06	6	250 ml	150 W (x6)
EAM-9203-06	6	500 ml	260 W (x6)
EAM-9204-06	6	1000 ml	350 W (x6)
EAM-9205-06	6	2000 ml	450 W (x6)