

ES-Series (Rotamantle)

Controlled Stirrer Mantles

Specifications:

Case Construction:	Steel, with chemical resistant coating
Thermal Insulation:	Moulded silica powder
Heating Element:	Nichrome wire
Mantle:	Woven glass fibre mat
Element Temperature:	Up to 450°C (Max)
Speed Range:	100 ~ 1000 rpm
Power Supply:	220V AC



MODEL NO.	CAPACITY	WATTS
ES-3005	50ml	60 W
ES-301	100ml	80 W
ES-302	250ml	150 W
ES-303	500ml	260 W
ES-304	1000ml	350 W
ES-305	2000ml	450 W
ES-306	3000ml	550 W
ES-307	5000ml	700 W
ES-308	10000ml	1200 W
ES-309	20000ml	1800 W