

# DIGIT-090

Pocket refractometer measuring high sugar concentrations

## Specifications:

**Scale:** (triple)

- 0-42% Brix
- 42-71% Brix
- 71-90% Brix

Adjustment ring for colour dispersion

Main glass prism fixed with a chemical resistant epoxy

Prism cover: plexi

Focusing eyepiece with hard plastic eyecup

Dimensions: 55 x 125 x 235mm.

Weight: 900g

Calibration: zero calibration with distilled water at 20°

**Precision**

0,2 % Brix

0,2 % Brix

0,2 % Brix

