



Scale 1:1

Order number 51343160
former order number -

Specifications
pH range 0 - 14
zero point (pH) 7.00 ± 0.25

temperature range 0 - 80 °C
temperature probe -

shaft material glass
membrane glass U
membrane resistance < 600 MΩ

reference system ARGENTHAL™, Ag⁺ trap
reference diaphragm ceramic
reference electrolyte 3 mol/L KCl

bridge electrolyte -

electrode head S7
cable length -
connector -

cell constant (cm⁻¹) -
cell type -

indicating metal -

Storage 3 mol/L KCl

Special features
micro electrode 3 mm
ARGENTHAL™ with Ag⁺ trap reference system

Choice of cables
This electrode has a grey S7 (InLab) connector. The cable selector guide offers a large choice of grey S7 cables. The grey S7 connectors are compatible with the red S7 connectors. The color is only design-related.

replaces the
InLab®423