

Cleaning and Disinfection

90cm wide Freestanding/built-under units

Performance specifications

Narrow-necked glassware: 37/E 414 TA

Pipettes: 96/E 408 TA

Test tubes: 160/E 103-139

Control system

MTN - Multitronic Novo Plus

Electronic controls with 8 programmes

Various process parameters can be modified

Display of programme sequence and faults

Temperature and countdown display

Features

Stainless-steel cabinet with 2 levels

Spray arms/Stainless steel

Docking system for upper baskets/mobile injector units

Quadruple filtration system

UP Qmax circulation pump 400 l/min

AP Drain pump

ET Electric door lock

PASI Programme interruption backup

SU Buzzer (acoustic signal)

SST Interface for process documentation

TA Drying unit/radial fan

WES Water softener

Depending on model, with/without:

ADP AD demineralised water pump

DK Steam condenser/aerosol injection

Heating

Electric: 6kW (230V) or 9kW(400V)

Dispensing technology

1 door dispenser for powder detergent

1 door dispenser for liquid products

1 dispenser pump 10/30 for liquid acidic agents

1 dispenser pump 60/30 for liquid detergent

Cleaning technology

Fresh water system

Thermal disinfection at 93°C with 10 min. holding time

Docking system for upper baskets and mobile injector units

Design

Frontloading unit with bottom-hinged door

Freestanding/built-under unit without lid

Slanted fascia panel

Outer casing:

AE Stainless-steel

H xW xD: 820 x 900 x 700 mm

Test certificates

VDE, GS, EMV/

Interference suppression

, DVGW, WRAS

G 7883CD Mielabor

