

# SONICATOR® 3000

The Sonicator® 3000 offers a powerful generator with a 600-Watt output, a 20kHz convertor, temperature control, auto-tuning, 9 level microprocessor-controlled memory, external cooling device interface, digital display, wrench set and quick-disconnect cable.

Part # S3000 is the Sonicator 3000® package that includes one each of the following:

- Generator
- Convertor
- 1/2" Horn
- 6' Cable
- Wrench Set
- Temperature probe

## SPECIFICATIONS

<b>Generator</b>	
Dimensions	8.5" wide x 19.0" deep x 10.5" high (21.6 cm x 48.3 cm x 26.7 cm)
Weight	17.6 lbs. (8.0 kg)
Maximum Power Output	600 Watts
Operating Frequency	20 kHz
Input	110VAC @ 10 Amps/220VAC @ 5 Amps
Total Processing Time	1 sec. to 100 hrs.
Adjustable Pulse On & Off Time	0.5 sec. to 1 hr.
<b>Convertor</b>	
Dimensions	2.5" diameter x 6.0" length (6.4 cm x 15.2 cm)
Weight	1.75 lbs. (0.79 kg)

