

# MICROPLATE HORN

The Microplate Horn comes complete with a required acoustic enclosure to protect users from the sound waves generated by this horn, as well as all the attachments required to control and replace the water in the horn reservoir.

Part # 431MP includes one each of the following:

- Microplate Horn
- Acoustic Enclosure
- Syringe
- Pinch Clamps
- Silicone Tubing

**Note:** The Microplate Horn is also available as part of a complete package, part # S3000MP which also includes the S3000 generator, convertor, cable and wrench set.



## SPECIFICATIONS

### Horn Reservoir

Dimensions	6.5" O.D. x 5.9" I.D. x 0.6" deep (16.5 cm x 15.0 cm x 1.5 cm)
------------	---

### Acoustic Enclosure

External Dimensions	10.0" wide x 10.0" deep x 17.0" high (25.4 cm x 25.4 cm x 43.2 cm)
---------------------	---