


## IN-LINE PROCESSING



## ORDERING INFORMATION

**800A** and **800B** are supplied complete with casing, horn,  $\frac{3}{32}$ " orifice plate and  $\frac{1}{8}$ " orifice plate.

### Replacement Orifice Plates

**Part #**

- 842**  $\frac{3}{32}$ " Orifice plate
- 843**  $\frac{1}{8}$ " Orifice plate
- 840** Blank orifice plate with pilot hole (extremely low flow rates)

**800C** is supplied as FloCell® casing only. Horn and matching flow tube must be specified and ordered separately.

**Part #**

- 821CT**  $\frac{3}{4}$ " Tapped processing horn
- 821C**  $\frac{3}{4}$ " Solid processing horn
- 822CT** 1" Tapped processing horn
- 822C** 1" Solid processing horn
- 845** Flow tube for use with 821C or 821CT
- 846** Flow tube for use with 822C or 822CT

**800D** is supplied as FloCell® casing with flow tube. The horn must be specified and ordered separately.

**Part #**

- 822DT** 1" Tapped horn
- 822D** 1" Solid horn

### Replacement Flow Tube

**Part #**

- 847** Flow tube

| Part #                 | 800A                      | 800B                | 800C                | 800D              |
|------------------------|---------------------------|---------------------|---------------------|-------------------|
| Maximum Flow Rating    | 0.67 lpm                  | 0.67 lpm            | 20 lpm              | 40 lpm            |
| Material of Body Tube  | Translucent Polycarbonate | Stainless Steel     | Stainless Steel     | Stainless Steel   |
| Static Pressure Rating | 50 psi (350 kPa)          | 100 psi (700 kPa)   | 100 psi (700 kPa)   | 100 psi (700 kPa) |
| Inlet/Outlet size      | $\frac{1}{4}$ " NPT       | $\frac{1}{4}$ " NPT | $\frac{1}{2}$ " NPT | 1" NPT            |

\*lpm = liters per minute