


# MISCELLANEOUS

## SONABOX™ ACOUSTIC ENCLOSURE

The Sonabox™ reduces noise down to safe, comfortable levels. All Sonicator® horns and attachments are compatible with the Sonabox™ except the 431MP which comes with a custom acoustic enclosure.



**Part # 432B**

Interior Dimensions: 9.25" Wide x 10.75" Deep x 19.75" High  
(23.5 cm x 27.3 cm x 50.2 cm)

## COOLING CELLS



**Part # 428**

Volume 0.5 to 10 ml

Water jacketed sonication vessel with conical well for use with Microtip probes.

**Part # 428M**

Volume 10 to 25 ml (as show above)

Water jacketed sonication vessel with flat-bottom well for use with 1/2" horns.

## DUAL HORN

The Dual Horn allows a single Sonicator® and convertor to process two samples simultaneously.



Suggested Sonicator® accessories for Environmental Testing (see EPA Method SW846-3550).

**Part # 351**

Dual Horn

**Part # 355**

3/4" Stepped Extenders (must order 2)

**Part # 356**

Dual Beaker Holder (optional)

**Part # 438**

Convertor Stand and Positioner (optional)

## ROSETT™

Conical shaped glass vessels with semicircular side arms to promote circulation and efficient cooling of sample. Available with jacket (J) for continuous coolant flow or without a jacket (call customer service for further information).

**Part # 425 & 425J**

Volume 5 to 15 ml

**Part # 427 & 427J**

Volume 100 to 450 ml