

NALGENE Mr. Frosty Freezing Container



NALGENE®

Mr. Frosty Freezing Container *Cooling Rate of -1° Per Minute*

Provides the critical, repeatable cooling rate required for cell cryopreservation and recovery. Requires only 100% isopropyl alcohol and mechanical

freezer. Containers with alcohol can be stored at room temperature, saving freezer space.

Specifications: Made of polycarbonate. White vial holder and blue closure made of high-density polyethylene. Labeled with graphic, step-by-step instructions. Foam insert. Holds up to 18 vials. Holder prevents vials from contacting alcohol--no contamination by wicking, no loss of labels or printing on vials.