

WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk Website: www.wolflabs.co.uk

Nichipet Ex

The Nichipet EX series of digital pipettes offers unsurpassed accuracy and light weight in a classical pipette profile. Eight models cover pipetting from 0.1µl up to 10ml. All models are completely autoclavable to ensure contamination-free pipetting.

- UV resistant body
- Unsurpassed accuracy
- Robust yet lightweight
- Fully autoclavable - minimises contamination risks
- Accurate digital volume setting
- Easy maintenance and calibration routines
- Thermally insulated from hand heat
- Filtered nozzles in 1000, 5000 and 10,000µl models prevent contamination
- Ceramic piston in 200, 1000, 5000 and 10,000µl models



Digital volume display
is highly visible and easily read by both right and left handed users



Volume setting via Thumb Knob
ensures rapid volume change (without problem of trapped gloves!!)



Volume locking lever
ensures no accidental changes to setting (what you set is what you get!!!)



NEW contour

the enlarged finger-rest provides increased comfort and extended fatigue free operation.

Tip Ejector Button

Convenient location eliminates the possibility of accidental tip ejection.

NEW UV resistant body

the NICHIPET-EX is now made in a material with enhanced UV resistance thereby making it the ideal pipette for use inside PCR clean cabinets.

Effortless Operation

Smooth, effortless movement of the piston reduces operator fatigue.

PTFE Seal / O-ring

New configuration ensures better air tightness, and better chemical resistance.

Full maintenance and calibration

service is available for all Nichiryo Liquid Handling products