

Optimax Microscopes

product datasheet

XT

Stereo zoom
microscope



Features and Specifications

- New ergonomic styling reduces strain during long or frequent periods of use
- High quality optics make observation of fine details simple
- Built-in base and incident light illumination
- Adjustable illumination intensity
- Built-in zoom range 0.7x to 4.2x
- Pillar support stand for easily adjustable height
- Adjustable interocular distance for comfortable viewing
- 10x widefield eyepieces
- Working distance typically 90mm
- Range of stands and bases available for different applications

Representing the very latest design in stereo microscopy. This ergonomic microscope includes high quality objectives whose clarity and brightness produce outstanding results. The XTL-101 includes zoom range encompasses 0.7x — 4.2x magnifications with a conveniently positioned zoom control knob. The integrated illumination in the base and incident illuminator maximises the types of specimens that can be observed.

Supplied to you by: