


Buffer Puffer™ self-recirculating system, A5



- ? Convenience of Gator
- ? Self recirculating buffer
- ? No uneven migration, band distortion or dissociation
- ? No external pumps, tubing or stirrer bars needed
- ? Ideal for long runs, RNA gels, multiple sample runs

Gel size 200mm(w)×250mm(l), buffer volume 2L. With buffer chamber, Super-Safe lid with attached power cord, gasketed UVT gel tray and 3 combs 1.5mm thick (16, 24, 36 wells).

Cat. No.	Each
306/7320/00	£471.73N

Accessories as Gator System, page 926

Product Code	Description	Pack Size
306732000	Gel system,self-recirculating,Buffer Puffer model A5	EACH