



● STR-1-LAB



● STR-2-LAB



● STR-3-LAB

Dimensions

Back : 400 x 420(h) mm
 Seat : 465 x 465 mm

Sitting Height :

STR-1-LAB : 605-855 mm
 STR-2-LAB : 530-730 mm
 STR-3-LAB : 460-590 mm



Technical Specifications

Chair Base :	Five star polished base with epoxy powder coated underside - Diametre 620 mm
HEPA Filter :	P2 complying with EN143:2000
Footring :	Chrome footring with polished ring and black epoxy powder coated underside - Diametre 460 mm Fully Adjustable
Castors/Glides :	Manufactured in accordance with DIN EN 12 528, ANSI BIFMA x 5.1 Load per castor 100N (DIN 51953)
Mechanism :	Independent seat and back tilt with back height adjustment
Loading Limit :	Gaslift maximum working limit 150 Kgs Castors/glides working limit 50 Kg each



Upholstery

Standard (Air and vapor permeability) or Medical Grade Vinyl (impermeable to air and water)
 Flame retardant to US MVSS 302 (inspected DIN 75 200)
 Flame retardant to BS 5852 Part 1