



## STANDARD REGULATOR & SELF SEALING VALVE

- Bull Nose Cylinder Fitting
- Brass Chrome Plated Body
- HP Gauge
- BS Gas Indexed Self Sealing Valve

R1601	Oxygen
R1611	MA-7
R1621	MA-4
R1601N	N <sub>2</sub> O, 11/16

R1601



## STANDARD REGULATOR & SELF SEALING VALVE

- Pin-Index Cylinder Fitting
- Brass Chrome Plated Body
- HP Gauge
- BS Gas Indexed Self Sealing Valve

R16015	Oxygen
R16115	MA-7
R16215	MA-4
R16015N	Nitrous Oxide
R16105/E	Entonox

R16015



## STANDARD REGULATOR PRE-SET

- 11/16 "BODOK" Cylinder Fitting
- Brass Chrome Plated Body
- Pre-Set 4Bar
- 1/8bsp male Outlet (non return)

R1605	Oxygen Nitrous Oxide
R1605L	Oxygen L/H
R1605L-N	N <sub>2</sub> O, L/H

R1605



## STANDARD REGULATOR PRE-SET

- 11/16 BULL NOSE Cylinder Fitting
- Brass Chrome Plated Body
- Pre-Set 4Bar
- 1/8bsp male Outlet (non return)

R1605-BN	Oxygen
R1605N-BN	Nitrous Oxide
R1605L-BN	Oxygen L/H
R1605LN-BN	N <sub>2</sub> O, L/H

R1605-BN



## STANDARD REGULATOR CO<sub>2</sub> ADJUSTABLE

- 14TPI Female Fitting
- Brass Chrome Plated Body
- HP Gauge
- 2nd Stage Regulator & 0 to 15psi Output Gauge
- Fine Adjustable Output Pressure
- Tubing Outlet

R1606	CO <sub>2</sub>
-------	-----------------

R1606



## STANDARD REGULATOR CO<sub>2</sub> ADJUSTABLE

- Pin-Index Cylinder Fitting
- Brass Chrome Plated Body
- HP Gauge
- 2nd Stage Regulator & 0 to 15psi Output Gauge
- Fine Adjustable Output Pressure
- Tubing Outlet

R16065	CO <sub>2</sub>
--------	-----------------

R16065