

WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk Website: www.wolflabs.co.uk

KF-75 ICE FLAKE MACHINE



Dimensions (mm) -

Width Depth Height

485 570 720

height with legs 860

Porkka KF 75		
Air cooled / Water cooled	Air	water
Production		
kg/24hr	75	78
Lb/24hr	165	172
Bin Capacity	25kg/55Lb	
Power (T=43C)-W	520	
Net Weight	54kg	
Gross Weight	62kg	
Volume When Packed	0.28 m 3	
Refrigerant used	R134a	
Standard Voltage	V230/50Hz/1	