

# WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: [sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)  
Website: [www.wolflabs.co.uk](http://www.wolflabs.co.uk)

## Compact X4 Automatic X-ray Film Processor



*The X4 processor is specifically set up by qualified staff to process the type of film emulsion used in a laboratory*

### **X-4 Processor Specification**

Cycle time – 120 seconds  
Film capacity – 70 films per hour  
Measured replenishment  
Comes complete with 25l replenishment tanks  
Microprocessor controlled  
Light tight feed-in tray

Complete system training is given after installation