

WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk Website: www.wolflabs.co.uk

Porkka KF 165 Ice Flake Machine



Porkka KF 165		
Air cooled / Water cooled	Air	water
Production		
kg/24hr	163	166
Lb/24hr	359	366
Bin Capacity	60kg/132Lb	
Power (T=43C)-W	650	
Net Weight	81kg	
Gross Weight	94kg	
Volume When Packed	0.64 m ³	
Refrigerant used	R134a	
Standard Voltage	V230/50Hz/1	

Dimensions (mm) -

Width Depth Height
955 590 780

height with legs - 920