

## INFECTION CONTROL SYSTEM TECHNICAL SPECIFICATION



<b>Product name:</b>	Infection Control System		
<b>Product code:</b>	739062	739064	739065

### COMPONENTS INCLUDED

Stainless steel tower housing unit	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
C3 Advance 16L vacuum autoclave	<input checked="" type="checkbox"/>	-	<input checked="" type="checkbox"/>
UltraClean thermal washer, disinfectant & dryer	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	-
Water purification system	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Water softener system	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

### DIMENSIONS

Width (m)	0.9
Depth (m)	0.7
Height (m)	1.5
Area required (m) (width x depth)	1 x 1.5
Max net weight (kg)	200

### INSTALLATION

Electricity supply Century autoclave (amp)	13
Electricity supply UltraClean (amp)	13
Electricity supply if 2 x printers supplied (amp)	2 x 5
Electricity supply if 2 x dataloggers supplied (amp)	1 x 5
Electricity supply water softener (amp)	5
Cold water supply pressure (bar)	0.5 – 3
Cold water supply temperature (°C)	15 – 25
Cold water supply hardness (mg/l)	0 – 210
Hot water supply pressure (bar)	0.5 – 3
Hot water supply temperature (°C)	50 – 60
Hot water supply hardness (mg/l)	0 – 210
Drainage (mm)	42mm trapped waste

### WARRANTY (Subject to product registration and based on individual components)

C3 Advance* - Parts (unlimited usage, excludes consumables)	1 Year
UltraClean - Parts (unlimited usage, excludes consumables)	1 Year

Prestige Medical will not be liable in the event that the purchaser has failed to adhere to instructions outlined in the user manual or if the equipment has been abused, interfered with, altered, repaired or serviced by any unauthorised party. This may result in the protection provided by the equipment being impaired. Consumer statutory rights are not affected. The Prestige Medical Ltd policy is one of continuous development and as such reserves the right to change the specification of models and items illustrated and described herein at any time. \*Terms and conditions apply, contact Prestige Medical for details.

Tel: +44 (0) 1254 682 622 Fax: +44 (0) 1254 682 606 Web: [www.prestigemedical.co.uk](http://www.prestigemedical.co.uk) Email: [info@prestigemedical.co.uk](mailto:info@prestigemedical.co.uk)