



Established for more than twenty years, Prestige Medical has been class leader in the development and manufacture of high technology equipment for sterilizing and disinfection.

Our equipment is fully approved and certified to the latest rigorous international standards, conforming to all the requirements of medical and legislative bodies.

With a reputation for working to world-class standards of performance and safety, our equipment has been sold to more than 80 countries - exceeding one million units delivered.

Prestige Medical continues to improve its products and systems for the increasingly demanding international market where our business continues to grow. Through our established technical strength and experience, as a full service support organisation, we help our customers develop and improve their infection control practices.

Our objective is consistently to provide the very best solutions in our sector for health practitioners across the world by ensuring that our specialist knowledge is applied in the most effective way.

PO Box 154,
Clarendon Road,
Blackburn, Lancashire,
BB1 9UG, England

Tel: +44 (0) 1254 682 622
Fax: +44 (0) 1254 682 606

www.prestigemedical.co.uk
sales@prestigemedical.co.uk

Registered in England
Reg No. 2826793
Vat No. 633707835

Registered Office:
Unit 1, First Avenue,
Maybrook Industrial Estate,
Minworth, Sutton Coldfield,
West Midlands, B76 1PA

Prestige Medical is a registered trademark.

D30 Washer Disinfector



D30 Washer Disinfector



key features

- Prestige Medical's new washer disinfector is designed for the cleaning and thermal disinfection of all types of instruments.
- 15 minute wash only cycle.
- 35 minute full wash and disinfect cycle.
- measured detergent dosing system - running costs only 6p per cycle.
- 8 dental handpiece internal wash system.
- built in detergent compartment and dispenser as standard.
- printer or data logger options.
- easy instrument loading and transfer system. holds up to 10 instrument cases.
- internal wash system for tattoo tubes.
- made entirely of stainless steel AISI316.
- easy programming with temperature and status displayed at any given moment.
- wide range of accessories.



providing the most effective solutions for instrument washing and disinfection...



technical specification

dimensions

width	500mm
depth	575mm
height	535mm
net weight	43kg
gross weight	46kg

connections

drain connection	dia 25mm
electric connection	230V/50Hz 400V/3+N/50Hz 415V/3+N/50Hz
power (Max)	2950w
pump power	150w
hydraulic connection cold water	3/4" gas
min/max pressure	1.0 - 5.0 bar
maximum capacity	24L/min at 2 bar pressure dynamic
min/max temperature	10° - 15°C
noise	52 dB (A)
min/max ambient temperature	+5°C/+40°C

standards and guidelines

designed and manufactured to comply with the following standards and guidelines:

- EN55011
- EN50082-1
- EN1441
- EN61010-1
- EN61010-2-45
- prEN15883-1

leading the way for innovation, reliability and customer service

furniture options




279543



729011



729012



739000



739003



739004



739005




739006




739007



739008



739010



739011