



## WHO

Primarily: Dentists

## WHAT

Prestige Medical UltracClean, Thermal Washer Disinfector with Drying

## WHY

1. Automate the washing of dental instruments
2. Comply with expected guidelines
3. Government grants available in their area specifically for washer disinfectors

## USPs

Offers **unique features** that **ensure continuing compliance** when pending guidelines are announced\*

Expected new guideline to recommend practices:

- replace manual washing with a thermal washer disinfector.
- integral drying of instruments
- In-built process verification

Also . . .

- Available as part of an infection control tower: a compact, free standing unit with Prestige Medical Autoclave

\* Money-back guarantee if the client can identify the UltraClean's non-compliance to any relevant new standards.

## TECHNICAL INFORMATION:

<b>Dimensions (mm):</b>	520 wide x 500 high x 490 deep
<b>Weight (kg):</b>	28 nominal
<b>Chamber :</b>	420 wide x 420 deep x 250 high
<b>dimensions (mm)</b>	
<b>Installation services:</b>	Hot water supply (if required) Cold water supply Drain Electricity Stable, secure worktop
<b>Consumption:</b>	Electricity 2.5kw Water: 6 litres/ cycle Detergent: 15ml cycle

## KEY BENEFITS:

- One of the **fastest cycles on the market**, with just 30 minutes wash, thermal disinfection and drying process time
- **Total compliance with the latest regulations** EN ISO15883 and NHS guidelines HTM2030 as well as the water regulations EN1717
- Quality guaranteed with a **1 year warranty**
- **Free-up staff** with an automated process that delivers consistent, effective results every time.
- **Easy to use** with a single programme and automatic detergent dispensing for optimum performance
- **Unique control system** that independently monitors all critical stages of the wash and disinfection process that proves compliance and provides alerts to faults.
- A choice of cold, RO (reverse osmosis) or hot water inlets to **speed up cycles**
- **Compact unit** that fits onto standard bench-tops and is simple to connect to domestic services.