

WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk

Website: www.wolflabs.co.uk

Lab Jacks



Manufactured to close tolerances from aluminium with stainless steel operating screw. The top and bottom plates are specially designed aluminium extrusion giving exceptional rigid and robust support without excessive weight. The plates are blue anodised finish. The labjacks are designed to be completely stable and free from rocking movements. Labjacks are supplied in blue.

All the labjacks with the exception of the 2041 and 2038 have a 10mm tapped hole in the right rear corner to take standard retort rods for which a locking nut is provided. Extension plates can be fitted to the 2036 & 2040 models.

| Model No | Platform Dimensions (mm) | Maximum Height (mm) | Closed Height (mm) | Gross Weight (kg) | Maximum Load Stationary (kg) | Maximum Operating Load (kg) |
|----------|--------------------------|---------------------|--------------------|-------------------|------------------------------|-----------------------------|
| 2041 | 50 x 40 | 78 | 36 | 0.16 | 20 | 10 |
| 2038 | 100 x 100 | 140 | 43 | 0.56 | 30 | 10 |
| 2035 | 150 x 150 | 275 | 60 | 1.4 | 60 | 25 |
| 2036 | 200 x 200 | 290 | 57 | 2.26 | 80 | 25 |
| 2040 | 200 x 230 | 290 | 57 | 2.42 | 80 | 25 |
| 2042 | 320 x 260 | 470 | 128 | 6.7 | 90 | 30 |

A range of non-slip mats are available. They are made of a polymeric compound.
Specifications

| Product Code | Size (mm) |
|--------------|-----------|
| 2047/1 | 50 x 40 |
| 2047/2 | 100 x 100 |
| 2047/4 | 150 x 150 |
| 2047/6 | 200 x 200 |
| 2047/7 | 200 x 230 |
| 2047/9 | 320 x 260 |

Also available is an extension plate for use with the 2036 & 2040 jacks. This will extend the platform dimensions to 300 x 300mm. Product code 2039.