



MOBILE FLAMMABLE CONTAINER

Steel container mounted on castors with lockable hinged lid and finished in a white epoxy coating. It has a removable steel insert with dividers for 8 x 2.5 liter bottles. The insert has leakproof seams and will retain spillages from all eight bottles. The container comes complete with Flammable Liquid symbol, index card holder and two keys.

A towing handle is available as an optional extra (Cat Ref E43.44).

Cat Ref	External Dimensions	Capacity
E43.40	350 x 670 x 400mm high	8 x 2.5 liter Euro/Winchester bottles



FLAMMABLE CABINET

Upright steel cabinet with lockable hinged door finished in a white epoxy coating. It has two levels with removable steel insert trays each with the capacity to hold 4 x 2.5 liter bottles. The insert trays have leakproof seams and will retain the spillage from one bottle. The cabinet comes complete with Flammable Liquid symbol, index card holder and two keys.

Cat Ref	External Dimensions	Capacity
E43.39	380 x 325 x 785mm high	8 x 2.5 liter Euro/Winchester bottles

Both of the above flammable liquid containers are designed to conform to the Highly Flammable & Liquefied Petroleum Gases Regulations 1972



SPECIMEN BUCKETS

White molded polythene buckets with press-on lid and carrying handle.

Cat Ref	Volume
E89.11 (without handle)	0.5 liter
E89.12	1.0 liter
E89.13	2.5 liter
E89.14	5.0 liter
E89.16	10.0 liter
E89.17	15.0 liter