

POSTLIP BLOTTING PAPER

Highly absorbent medium fiber-free blotting paper. Supplied in boxes of 100.

Cat Ref	Size
E1	89 x 140mm



FORDS GOLD MEDAL BLOTTING PAPER

Highly absorbent fiber-free blotting paper. Supplied in packs of 25 sheets.

Cat Ref	Size
E1.01	445 x 570mm



WHATMAN NO 1 FILTER PAPER CIRCLES

Quality Whatman No 1 filter paper circles. Supplied in boxes of 100.

Cat Ref	Diameter	Cat Ref	Diameter
E10/F1	7.0cm	E10/F6	18.5cm
E10/F2	9.0cm	E10/F7	24.0cm
E10/F3	11.0cm	E10/F8	27.0cm
E10/F4	12.5cm	E10/F9	32.0cm
E10/F5	15.0cm	Other sizes available upon request.	



WHATMAN FILTER PAPER

Quality Whatman No 1 filter paper. Supplied in boxes of 100 sheets.

Cat Ref	Size
E10/F10	46 x 57cm
E10/F10.1	58 x 68cm

Other sizes available upon request.

WHATMAN BENCHKOTE AND BENCHKOTE PLUS

Whatman Benchkote is a high quality white absorbent paper coated on one side with an impermeable layer of polypropylene. Used as a disposable bench covering.

Cat Ref	Size	Pack
E10.97/4200	46 x 57cm	50 sheets
E10.97/4201	46cm x 50m	Reel
E10.97/4202	92cm x 50m	Reel

Whatman Benchkote Plus is a high quality white paper with high absorbency.

Cat Ref	Size	Pack
E10.97/6150	50 x 60cm	50 sheets
E10.97/6160	60cm x 50m	Reel

